

Antibodies Product Datasheet

Product Name:	Cytokeratin 19 Monoclonal Antibody
Catalog Number:	N0185
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Cytokeratin 19 expressed in E. Coli.
Concentration:	
Gene Name:	KRT19
Protein Name:	Keratin type I cytoskeletal 19
Other Name:	KRT19; Keratin; type I cytoskeletal 19; Cytokeratin-19; CK-19; Keratin-19; K19
Human Swiss Prot No.:	P08727
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Cytokeratin 19 Monoclonal Antibody detects endogenous levels of Cytokeratin 19 protein.
Storage and Stability:	-20°C/1 year