

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Luciferase mouse Monoclonal Antibody(7C9)
<b>Catalog Number:</b>	N1639
<b>Reactivity:</b>	Firefly
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:1000-2000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthetic Peptide of Luciferase
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	
<b>Protein Name:</b>	Luciferase
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	P08659
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Luciferase protein detects endogenous levels of Luciferase
<b>Storage and Stability:</b>	-20°C/1 year