
Antibodies Product Datasheet

Product Name:	Smad3 mouse Monoclonal Antibody(1E8)
Catalog Number:	N1633
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Smad3 at AA range of 350-430
Concentration:	1 mg/ml
Gene Name:	SMAD3
Protein Name:	SMAD3
Other Name:	SMAD3
Human Swiss Prot No.:	P84022
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Smad3 protein detects endogenous levels of SMAD3
Storage and Stability:	-20°C/1 year