

Antibodies Product Datasheet

Product Name:	CAMK1D Antibody
Catalog Number:	C49952
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CAMK1D
Concentration:	1 mg/ml
Gene Name:	CAMK1D CAMKID
Protein Name:	Calcium/calmodulin-dependent protein kinase type 1D (EC 2.7.11.17) (CaM kinase I delta) (CaM kinase ID) (CaM-KI delta) (CaMKI delta) (CaMKID) (CaMKI-like protein kinase) (CKLiK)
Other Name:	
Human Swiss Prot No.:	Q8IU85
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CAMK1D at Human, Mouse
Storage and Stability:	Store at -20° C/1 year