

Antibodies Product Datasheet

Product Name:	IkappaB-beta Antibody
Catalog Number:	R12-2939
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	IkappaB-beta Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IkappaB-beta.
Concentration:	1 mg/ml
Gene Name:	NFKB1B
Other Name:	I-kappa-B-beta; IKB-B; IKBB; IkappaBbeta; NF-kappa-BIB; NF-kappaB inhibitor beta; NFKB1B; TR-interacting protein 9; TRIP9; Thyroid receptor interacting protein 9
Human Swiss Prot No.:	Q15653
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	IkappaB-beta Antibody antibody detects endogenous levels of IkappaB-beta.
Storage and Stability:	Stable at -20°C for at least 1 year.