

## Antibodies Product Datasheet

<b>Product Name:</b>	FEN-1 mouse mAb
<b>Catalog Number:</b>	N0978
<b>Reactivity:</b>	Human, Mouse, Rat, Monkey, Hamster
<b>Applications:</b>	WB, ICC
<b>Dilutions:</b>	wb dilution 1:1000 icc dilution 1:400
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Purified recombinant human FEN-1 protein fragments expressed in E.coli.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	fen1
<b>Protein Name:</b>	
<b>Other Name:</b>	DNase IV; FEN-1; FEN1; FEN1_HUMAN; Flap endonuclease 1; Flap structure specific endonuclease 1; Flap structure-specific endonuclease 1; hFEN-1; hFEN1; Maturation factor 1; MF1; Rad2.
<b>Human Swiss Prot No.:</b>	P39748
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of FEN-1 and does not cross-react with related proteins.
<b>Storage and Stability:</b>	-20°C/1 year