
Antibodies Product Datasheet

Product Name:	OR6K2 Antibody
Catalog Number:	R12-3253
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	OR6K2 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human OR6K2.
Concentration:	1 mg/ml
Gene Name:	OR6K2
Other Name:	Olfactory receptor 6K2; Olfactory receptor OR1-17; OR6K2
Human Swiss Prot No.:	Q8NGY2
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	OR6K2 Antibodyantibody detects endogenous levels of OR6K2.
Storage and Stability:	Stable at -20°C for at least 1 year.