
Antibodies Product Datasheet

Product Name:	Cleaved-CD97 β (S531) Antibody
Catalog Number:	L0192
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD97 beta. AA range:512-561
Concentration:	1 mg/ml
Gene Name:	CD97
Protein Name:	CD97 antigen
Other Name:	CD97; CD97 antigen; Leukocyte antigen CD97; CD antigen CD97
Human Swiss Prot No.:	P48960
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-CD97 β (S531) Polyclonal Antibody detects endogenous levels of fragment of activated CD97 β protein resulting from cleavage adjacent to S531.
Storage and Stability:	-20°C/1 year