

Antibodies Product Datasheet

Product Name:	hnRNP F Antibody
Catalog Number:	C10700
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human hnRNP F. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	HNRNPF
Protein Name:	Heterogeneous nuclear ribonucleoprotein F
Other Name:	HNRNPF; HNRPF; Heterogeneous nuclear ribonucleoprotein F; hnRNP F; Nucleolin-like protein mcs94-1
Human Swiss Prot No.:	P52597
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	hnRNP F Polyclonal Antibody detects endogenous levels of hnRNP F protein.
Storage and Stability:	-20°C/1 year