

Antibodies Product Datasheet

Product Name:	Jagged1 Antibody
Catalog Number:	C30483
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human JAG1. AA range:981-1030
Concentration:	1 mg/ml
Gene Name:	JAG1
Protein Name:	Protein jagged-1
Other Name:	JAG1; JAGL1; Protein jagged-1; Jagged1; hJ1; CD339
Human Swiss Prot No.:	P78504
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Jagged1 Polyclonal Antibody detects endogenous levels of Jagged1 protein.
Storage and Stability:	-20°C/1 year