
Antibodies Product Datasheet

Product Name:	GSK3 β mouse Monoclonal Antibody(4D2)
Catalog Number:	N1617
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of GSK3 β at AA range of 1-80
Concentration:	1 mg/ml
Gene Name:	GSK3B
Protein Name:	GSK3B
Other Name:	GSK3B
Human Swiss Prot No.:	P49841
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GSK3 β protein detects endogenous levels of GSK3B
Storage and Stability:	-20°C/1 year