

Antibodies Product Datasheet

Product Name:	CSB Antibody
Catalog Number:	C10869
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ERCC6. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	ERCC6
Protein Name:	DNA excision repair protein ERCC-6
Other Name:	ERCC6; CSB; DNA excision repair protein ERCC-6; ATP-dependent helicase ERCC6; Cockayne syndrome protein CSB
Human Swiss Prot No.:	Q03468
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CSB Polyclonal Antibody detects endogenous levels of CSB protein.
Storage and Stability:	-20°C/1 year