

## Antibodies Product Datasheet

<b>Product Name:</b>	p52 S6 kinase Antibody
<b>Catalog Number:</b>	C11546
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human RPS6KC1. AA range:231-280
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	RPS6KC1
<b>Protein Name:</b>	Ribosomal protein S6 kinase delta-1
<b>Other Name:</b>	RPS6KC1; RPK118; Ribosomal protein S6 kinase delta-1; S6K-delta-1; 52 kDa ribosomal protein S6 kinase; Ribosomal S6 kinase-like protein with two PSK domains 118 kDa protein; SPHK1-binding protein
<b>Human Swiss Prot No.:</b>	Q96S38
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	p52 S6 kinase Polyclonal Antibody detects endogenous levels of p52 S6 kinase protein.
<b>Storage and Stability:</b>	-20°C/1 year