

Antibodies Product Datasheet

Product Name:	SVCT2 Antibody
Catalog Number:	C49975
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SVCT2
Concentration:	1 mg/ml
Gene Name:	SLC23A2 KIAA0238 NBTL1 SLC23A1 SVCT2 YSPL2
Protein Name:	Solute carrier family 23 member 2 (Na ⁺)/L-ascorbic acid transporter 2) (Nucleobase transporter-like 1 protein) (Sodium-dependent vitamin C transporter 2) (hSVCT2) (Yolk sac permease-like molecule 2)
Other Name:	
Human Swiss Prot No.:	Q9UGH3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SVCT2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year