

Antibodies Product Datasheet

Product Name:	PNMA2 Antibody
Catalog Number:	C42364
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human PNMA2
Concentration:	1 mg/ml
Gene Name:	PNMA2 KIAA0883 MA2
Protein Name:	Paraneoplastic antigen Ma2 (40 kDa neuronal protein) (Onconeural antigen Ma2) (Paraneoplastic neuronal antigen MM2)
Other Name:	
Human Swiss Prot No.:	Q9UL42
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human PNMA2
Storage and Stability:	-20°C/1 year