

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Cdk1 Monoclonal Antibody
<b>Catalog Number:</b>	N0143
<b>Reactivity:</b>	Human, Mouse, Rat, Monkey
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant fragment of human Cdk1 expressed in E. Coli.
<b>Concentration:</b>	
<b>Gene Name:</b>	CDK1
<b>Protein Name:</b>	Cell division protein kinase 1
<b>Other Name:</b>	CDK1; CDC2; CDC28A; CDKN1; P34CDC2; Cyclin-dependent kinase 1; CDK1; Cell division control protein 2 homolog; Cell division protein kinase 1; p34 protein kinase
<b>Human Swiss Prot No.:</b>	P06493
<b>Formulation:</b>	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
<b>Specificity:</b>	Cdk1 Monoclonal Antibody detects endogenous levels of Cdk1 protein.
<b>Storage and Stability:</b>	-20°C/1 year