

Antibodies Product Datasheet

Product Name:	CNDP1 Antibody
Catalog Number:	C43182
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 330-410
Concentration:	1 mg/ml
Gene Name:	CNDP1 CN1 CPGL2 UNQ1915/PRO4380
Protein Name:	Beta-Ala-His dipeptidase (EC 3.4.13.20) (CNDP dipeptidase 1) (Carnosine dipeptidase 1) (Glutamate carboxypeptidase-like protein 2) (Serum carnosinase)
Other Name:	
Human Swiss Prot No.:	Q96KN2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CNDP1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year