
Antibodies Product Datasheet

Product Name:	KIT (Phospho-Tyr936) Antibody
Catalog Number:	P12-1132
Reactivity:	Human, Mouse
Applications:	WB, IHC, ELISA
Dilutions:	KIT (Phospho-Tyr936) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human KIT around the phosphorylation site of Tyr936.
Concentration:	1 mg/ml
Gene Name:	KIT
Other Name:	SCFR; Proto-oncogene tyrosine-protein kinase Kit; c-kit; CD117 antigen
Human Swiss Prot No.:	P10721
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	KIT (Phospho-Tyr936) Antibody detects endogenous levels of KIT only when phosphorylated at Tyr936.
Storage and Stability:	Stable at -20°C for at least 1 year.