

Antibodies Product Datasheet

Product Name:	MMP-19 Antibody
Catalog Number:	C0269
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MMP-19. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	MMP19
Protein Name:	Matrix metalloproteinase-19
Other Name:	MMP19; MMP18; RASI; Matrix metalloproteinase-19; MMP-19; Matrix metalloproteinase RASI; Matrix metalloproteinase-18; MMP-18
Human Swiss Prot No.:	Q99542
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MMP-19 Polyclonal Antibody detects endogenous levels of MMP-19 protein.
Storage and Stability:	-20°C/1 year