

Antibodies Product Datasheet

Product Name:	POK2 Antibody
Catalog Number:	C45604
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human POK2
Concentration:	1 mg/ml
Gene Name:	ERVK6
Protein Name:	POK2
Other Name:	HERV-K_7p22.1 provirus ancestral Pol protein (HERV-K(C7) Pol protein) (HERV-K(HML-2.HOM) Pol protein) (HERV-K108 Pol protein) [Includes: Reverse transcriptase (RT) (EC 2.7.7.49); Ribonuclease H (RNase H) (EC 3.1.26.4); Integrase (IN)]
Human Swiss Prot No.:	Q9BXR3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.241% sodium azide.
Specificity:	This antibody detects endogenous levels of POK2 at Human
Storage and Stability:	-20°C/1 year