

Antibodies Product Datasheet

Product Name:	DHRS2 Antibody
Catalog Number:	C15322
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DHRS2. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	DHRS2
Protein Name:	Dehydrogenase/reductase SDR family member 2
Other Name:	DHRS2; Dehydrogenase/reductase SDR family member 2; Dicarbonyl reductase HEP27; Protein D
Human Swiss Prot No.:	Q13268
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DHRS2 Polyclonal Antibody detects endogenous levels of DHRS2 protein.
Storage and Stability:	-20°C/1 year