

Antibodies Product Datasheet

Product Name:	EXOC1 Antibody
Catalog Number:	C48168
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EXOC1
Concentration:	1 mg/ml
Gene Name:	EXOC1 SEC3 SEC3L1 BM-012
Protein Name:	Exocyst complex component 1 (Exocyst complex component Sec3)
Other Name:	
Human Swiss Prot No.:	Q9NV70
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EXOC1 at Human
Storage and Stability:	Store at -20° C/1 year