

Antibodies Product Datasheet

Product Name:	TAL1 Antibody
Catalog Number:	B8431
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TAL-1. AA range:96-145
Concentration:	1 mg/ml
Gene Name:	TAL1
Protein Name:	T-cell acute lymphocytic leukemia protein 1
Other Name:	TAL1; BHLHA17; SCL; TCL5; T-cell acute lymphocytic leukemia protein 1; TAL-1; Class A basic helix-loop-helix protein 17; bHLHa17; Stem cell protein; T-cell leukemia/lymphoma protein 5
Human Swiss Prot No.:	P17542
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TAL1 Polyclonal Antibody detects endogenous levels of TAL1 protein.
Storage and Stability:	-20°C/1 year