
Antibodies Product Datasheet

Product Name:	TNR16 Antibody
Catalog Number:	R12-3640
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	TNR16 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TNR16.
Concentration:	1 mg/ml
Gene Name:	NGFR
Other Name:	Gp80-LNGFR; low affinity neurotrophin receptor p75NTR; low-affinity nerve growth factor receptor; NGF receptor; p75 ICD; TNFRSF16; TNR16; tumor necrosis factor receptor superfamily member 16
Human Swiss Prot No.:	P08138
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	TNR16 Antibodyantibody detects endogenous levels of TNR16.
Storage and Stability:	Stable at -20°C for at least 1 year.