

Antibodies Product Datasheet

Product Name:	ELOB Antibody
Catalog Number:	C43607
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	TCEB2
Protein Name:	Transcription elongation factor B polypeptide 2 (Elongin 18 kDa subunit) (Elongin-B) (EloB) (RNA polymerase II transcription factor SIII subunit B) (SIII p18)
Other Name:	
Human Swiss Prot No.:	Q15370
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ELOB Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year