

Antibodies Product Datasheet

Product Name:	Caspase-9 Antibody
Catalog Number:	B0832
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Caspase 9. AA range:119-168
Concentration:	1 mg/ml
Gene Name:	CASP9
Protein Name:	Caspase9
Other Name:	CASP9; MCH6; Caspase-9; CASP-9; Apoptotic protease Mch-6; Apoptotic protease-activating factor 3; APAF-3; ICE-like apoptotic protease 6; ICE-LAP6
Human Swiss Prot No.:	P55211
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Caspase-9 Polyclonal Antibody detects endogenous levels of Caspase-9 protein.
Storage and Stability:	-20°C/1 year