

Antibodies Product Datasheet

Product Name:	syntenin Antibody
Catalog Number:	C40075
Reactivity:	Human, Rat, Dog, Pig
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human syntenin
Concentration:	1 mg/ml
Gene Name:	SDCBP MDA9 SYCL
Protein Name:	syntenin
Other Name:	Syntenin-1 (Melanoma differentiation-associated protein 9) (MDA-9) (Pro-TGF- α cytoplasmic domain-interacting protein 18) (TACIP18) (Scaffold protein Pbp1) (Syndecan-binding protein 1)
Human Swiss Prot No.:	O00560
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human,Rat,Dog,Pig syntenin
Storage and Stability:	-20°C/1 year