
Antibodies Product Datasheet

Product Name:	KLF2 Antibody
Catalog Number:	C44422
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KLF2
Concentration:	1 mg/ml
Gene Name:	KLF2 LKLF
Protein Name:	KLF2
Other Name:	Krueppel-like factor 2 (Lung krueppel-like factor)
Human Swiss Prot No.:	Q9Y5W3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.70% sodium azide.
Specificity:	This antibody detects endogenous levels of KLF2 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year