

Antibodies Product Datasheet

Product Name:	MEK-2 Antibody
Catalog Number:	B7149
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), IP, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAP2K2. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	MAP2K2
Protein Name:	Dual specificity mitogen-activated protein kinase kinase 2
Other Name:	MAP2K2; MEK2; MKK2; PRKMK2; Dual specificity mitogen-activated protein kinase kinase 2; MAP kinase kinase 2; MAPKK 2; ERK activator kinase 2; MAPK/ERK kinase 2; MEK 2
Human Swiss Prot No.:	P36507
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MEK-2 Polyclonal Antibody detects endogenous levels of MEK-2 protein.
Storage and Stability:	-20°C/1 year