
Antibodies Product Datasheet

Product Name:	Scribble (Phospho-Ser1220) Antibody
Catalog Number:	A9619
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Scribble (Ser1220)
Concentration:	1 mg/ml
Gene Name:	SCRIB CRIB1 KIAA0147 LAP4 SCRB1 VARTUL
Protein Name:	Scribble (Ser1220)
Other Name:	Protein scribble homolog (Scribble) (hScrib) (Protein LAP4)
Human Swiss Prot No.:	Q14160
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Scribble (phospho-Ser1220)
Storage and Stability:	-20°C/1 year