
Antibodies Product Datasheet

| | |
|-------------------------------|-------------------------------------------------------------------------------------------------|
| Product Name: | BIOTIN*Polyclonal Rabbit Anti Dog IgG(h+l) |
| Catalog Number: | SA0284 |
| Reactivity: | DogIgG(h+l) |
| Conjugate: | Biotin |
| Applications: | WB,IHC-p,IF/ICC,ELISA |
| Dilutions: | WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000 |
| Modification: | Not Modification Specific |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| Concentration: | 4-5mg/mL |
| Other Name: | BIOTIN*Polyclonal Rabbit Anti-Dog IgG(h+l) |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage and Stability: | -20°C/1 year |