

Antibodies Product Datasheet

Product Name:	Histone H2B Monoclonal Antibody(Mix)
Catalog Number:	N1427
Reactivity:	Human,Mouse,Rat
Applications:	WB,IF/ICC
Dilutions:	WB: 1:1000-3000 IF 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Histone H2B
Concentration:	
Gene Name:	HIST1H2BB
Protein Name:	Histone H2B type 1-B
Other Name:	HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B.f; H2B/f
Human Swiss Prot No.:	P33778
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Histone H2B protein.
Storage and Stability:	-20°C/1 year