

Antibodies Product Datasheet

Product Name:	p130 Cas (Phospho-Tyr128) Antibody
Catalog Number:	P12-1038
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	p130 Cas (Phospho-Tyr128) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human p130 Cas around the phosphorylation site of Tyr128.
Concentration:	1 mg/ml
Gene Name:	BCAR1
Other Name:	Breast cancer anti-estrogen resistance protein 1; CRK-associated substrate; Cas scaffolding protein family member 1; p130cas; BCAR1; CAS; CASS1; CRKAS
Human Swiss Prot No.:	P56945
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	p130 Cas (Phospho-Tyr128) antibody detects endogenous levels of p130 Cas only when phosphorylated at Tyr128.
Storage and Stability:	Stable at -20°C for at least 1 year.