

Antibodies Product Datasheet

Product Name:	p18 Antibody
Catalog Number:	C0288
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human p18 INK. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	CDKN2C
Protein Name:	Cyclin-dependent kinase 4 inhibitor C
Other Name:	CDKN2C; CDKN6; Cyclin-dependent kinase 4 inhibitor C; Cyclin-dependent kinase 6 inhibitor; p18-INK4c; p18-INK6
Human Swiss Prot No.:	P42773
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	p18 Polyclonal Antibody detects endogenous levels of p18 protein.
Storage and Stability:	-20°C/1 year