

## Antibodies Product Datasheet

<b>Product Name:</b>	MKP-1 (Phospho Ser323) Antibody
<b>Catalog Number:</b>	C51138
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC, WB
<b>Dilutions:</b>	WB 1:500-2000 IHC 1:50-200
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Rabbit, polyclonal
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human MKP-1 (Phospho Ser323)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	DUSP1 CL100 MKP1 PTPN10 VH1
<b>Protein Name:</b>	Dual specificity protein phosphatase 1 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase hVH1) (Mitogen-activated protein kinase phosphatase 1) (MAP kinase phosphatase 1) (MKP-1) (Protein-tyrosine phosphatase CL100)
<b>Other Name:</b>	Dual specificity protein phosphatase 1 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase hVH1) (Mitogen-activated protein kinase phosphatase 1) (MAP kinase phosphatase 1) (MKP-1) (Protein-tyrosine phosphatase CL100)
<b>Human Swiss Prot No.:</b>	P28562
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of MKP-1 (Phospho Ser323) Rabbit pAb at Human, Mouse, Rat
<b>Storage and Stability:</b>	Store at -20° C/1 year