
Antibodies Product Datasheet

Product Name:	Cytokeratin 10 mouse mAb(PT1859)
Catalog Number:	N1833
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	WB 500-2000 IHC-p 1:100-500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG1, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Cytokeratin 10
Concentration:	
Gene Name:	KRT10 KPP
Protein Name:	Cytokeratin 10
Other Name:	Keratin, type I cytoskeletal 10 (Cytokeratin-10) (CK-10) (Keratin-10) (K10)
Human Swiss Prot No.:	P13645
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Cytokeratin 10. Heat-induced epitope retrieval (HIER) Citrate buffer of pH6.0 was highly recommended as antigen repair method in paraffin section
Storage and Stability:	-20°C/1 year