

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | FoxE3 Antibody |
| Catalog Number: | C10979 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, IHC-p, IF(paraffin section), ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human FOXE3. AA range:81-130 |
| Concentration: | 1 mg/ml |
| Gene Name: | FOXE3 |
| Protein Name: | Forkhead box protein E3 |
| Other Name: | FOXE3; FKHL12; FREAC8; Forkhead box protein E3; Forkhead-related protein FKHL12; Forkhead-related transcription factor 8; FREAC-8 |
| Human Swiss Prot No.: | Q13461 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | FoxE3 Polyclonal Antibody detects endogenous levels of FoxE3 protein. |
| Storage and Stability: | -20°C/1 year |