

Antibodies Product Datasheet

Product Name:	SDCCAG10 Antibody
Catalog Number:	C50855
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SDCCAG10
Concentration:	1 mg/ml
Gene Name:	CWC27 SDCCAG10 UNQ438/PRO871
Protein Name:	Peptidyl-prolyl cis-trans isomerase CWC27 homolog (PPIase CWC27) (EC 5.2.1.8) (Antigen NY-CO-10) (Serologically defined colon cancer antigen 10)
Other Name:	
Human Swiss Prot No.:	Q6UX04
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SDCCAG10 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year