

Antibodies Product Datasheet

Product Name:	Zic1/2/3/4/5 Antibody
Catalog Number:	C10230
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZIC1/2/3/4/5. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	ZIC1
Protein Name:	Zinc finger protein ZIC 1
Other Name:	ZIC1; ZIC; ZNF201; Zinc finger protein ZIC 1; Zinc finger protein 201; Zinc finger protein of the cerebellum 1; ZIC2; Zinc finger protein ZIC 2; Zinc finger protein of the cerebellum 2; ZIC3; ZNF203; Zinc finger protein ZIC 3; Zinc finger p
Human Swiss Prot No.:	Q15915/O95409/O60481/Q8N9L1/Q96T25
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Zic1/2/3/4/5 Polyclonal Antibody detects endogenous levels of Zic1/2/3/4/5 protein.
Storage and Stability:	-20°C/1 year