

Antibodies Product Datasheet

Product Name:	BCL2 (Phospho-Ser70) Antibody
Catalog Number:	P12-1079
Reactivity:	M, R
Applications:	WB, ELISA
Dilutions:	BCL2 (Phospho-Ser70) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human BCL2 around the phosphorylation site of Ser70.
Concentration:	1 mg/ml
Gene Name:	BCL2
Other Name:	Apoptosis regulator Bcl-2; B-cell CLL/lymphoma 2; BCL2
Human Swiss Prot No.:	P10415
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	BCL2 (Phospho-Ser70) Antibody detects endogenous levels of BCL2 only when phosphorylated at Ser70.
Storage and Stability:	Stable at -20°C for at least 1 year.