

Antibodies Product Datasheet

Product Name:	RBM1A1 Antibody
Catalog Number:	C50051
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human RBM1A1
Concentration:	1 mg/ml
Gene Name:	RBM1A1 RBM1 RBM2 YRRM1 YRRM2
Protein Name:	RNA-binding motif protein, Y chromosome, family 1 member A1 (RNA-binding motif protein 1) (RNA-binding motif protein 2) (Y chromosome RNA recognition motif 1) (hRBM1A)
Other Name:	
Human Swiss Prot No.:	P0DJJ3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RBM1A1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year