

Antibodies Product Datasheet

Product Name:	HSP27 Monoclonal Antibody(Mix-mA)
Catalog Number:	N1524
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB: 1:1000-1:3000 IHC: 1:50-1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of HSP27
Concentration:	1 mg/ml
Gene Name:	HSPB1
Protein Name:	Heat shock protein beta-1 (HspB1) (28 kDa heat shock protein) (Estrogen-regulated 24 kDa protein) (Heat shock 27 kDa protein) (HSP 27) (Stress-responsive protein 27) (SRP27)
Other Name:	Heat shock protein beta-1 (HspB1) (28 kDa heat shock protein) (Estrogen-regulated 24 kDa protein) (Heat shock 27 kDa protein) (HSP 27) (Stress-responsive protein 27) (SRP27)
Human Swiss Prot No.:	P04792
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Antibody can detects endogenous HSP27 protein.
Storage and Stability:	-20°C/1 year