

Antibodies Product Datasheet

Product Name:	MYCBP2 Antibody
Catalog Number:	C49790
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MYCBP2
Concentration:	1 mg/ml
Gene Name:	MYCBP2 KIAA0916 PAM
Protein Name:	Probable E3 ubiquitin-protein ligase MYCBP2 (EC 6.3.2.-) (Myc-binding protein 2) (Pam/highwire/rpm-1 protein) (Protein associated with Myc)
Other Name:	
Human Swiss Prot No.:	O75592
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MYCBP2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year