## **Assay Biotechnology Company**

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

## **Antibodies Product Datasheet**

**Product Name:** Sp1 (Phospho-Thr739) Antibody

Catalog Number: A7218

Reactivity: Human, Mouse, Rat, Monkey

**Applications:** WB, IHC-p, IF(paraffin section), ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000. Not yet tested in other applications.

Modification: Phospho Specific

**Source:** Polyclonal Rabbit

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Immunogen:** The antiserum was produced against synthesized peptide derived from

human SP1 around the phosphorylation site of Thr739. AA range:706-

755

Concentration: 1 mg/ml

Gene Name: SP1

**Protein Name:** Transcription factor Sp1

Other Name: SP1; TSFP1; Transcription factor Sp1

Human Swiss Prot No.: P08047

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

**Specificity:** Phospho-Sp1 (T739) Polyclonal Antibody detects endogenous levels of

Sp1 protein only when phosphorylated at T739.

Storage and Stability: -20°C/1 year