
Antibodies Product Datasheet

Product Name:	Bmx Monoclonal Antibody
Catalog Number:	N0075
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Bmx expressed in E. Coli.
Concentration:	
Gene Name:	BMX
Protein Name:	Cytoplasmic tyrosine-protein kinase BMX
Other Name:	BMX; Cytoplasmic tyrosine-protein kinase BMX; Bone marrow tyrosine kinase gene in chromosome X protein; Epithelial and endothelial tyrosine kinase; ETK; NTK38
Human Swiss Prot No.:	P51813
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Bmx Monoclonal Antibody detects endogenous levels of Bmx protein.
Storage and Stability:	-20°C/1 year