

Antibodies Product Datasheet

Product Name:	FRS2- α (Phospho-Tyr196) Antibody
Catalog Number:	A9642
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human FRS2- α (Tyr196)
Concentration:	1 mg/ml
Gene Name:	FRS2
Protein Name:	FRS2- α (Tyr196)
Other Name:	Fibroblast growth factor receptor substrate 2 (FGFR substrate 2) (FGFR-signaling adaptor SNT) (Suc1-associated neurotrophic factor target 1) (SNT-1)
Human Swiss Prot No.:	Q8WU20
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse FRS2- α (phospho-Tyr196)
Storage and Stability:	-20°C/1 year