

Antibodies Product Datasheet

Product Name:	Stathmin(PT2148) mouse mAb
Catalog Number:	N1873
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:100-500, WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG2b, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Stathmin
Concentration:	
Gene Name:	STMN1 C1orf215 LAP18 OP18
Protein Name:	Stathmin (Leukemia-associated phosphoprotein p18) (Metablastin) (Oncoprotein 18) (Op18) (Phosphoprotein p19) (pp19) (Prosolin) (Protein Pr22) (pp17)
Other Name:	
Human Swiss Prot No.:	P16949
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Stathmin (Leukemia-associated phosphoprotein p18) (Metablastin) (Oncoprotein 18) (Op18) (Phosphoprotein p19) (pp19) (Prosolin) (Protein Pr22) (pp17). Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH9.0
Storage and Stability:	-20°C/1 year