

Antibodies Product Datasheet

Product Name:	ZAC1 Antibody
Catalog Number:	C20728
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PLAGL1. AA range:311-360
Concentration:	1 mg/ml
Gene Name:	PLAGL1
Protein Name:	Zinc finger protein PLAGL1
Other Name:	PLAGL1; LOT1; ZAC; Zinc finger protein PLAGL1; Lost on transformation 1; LOT-1; Pleiomorphic adenoma-like protein 1; Tumor suppressor ZAC
Human Swiss Prot No.:	Q9UM63
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ZAC1 Polyclonal Antibody detects endogenous levels of ZAC1 protein.
Storage and Stability:	-20°C/1 year