

## Antibodies Product Datasheet

<b>Product Name:</b>	IL-17F Antibody
<b>Catalog Number:</b>	C30629
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	IHC-p 1:50-200, ELISA 1:10000-20000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthetic peptide from human protein at AA range: 81-130
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	IL17F IL24
<b>Protein Name:</b>	Interleukin-17F (IL-17F) (Cytokine ML-1) (Interleukin-24) (IL-24)
<b>Other Name:</b>	Interleukin-17F (IL-17F) (Cytokine ML-1) (Interleukin-24) (IL-24)
<b>Human Swiss Prot No.:</b>	Q96PD4
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	The antibody detects endogenous IL-17F
<b>Storage and Stability:</b>	-20°C/1 year