

## Antibodies Product Datasheet

<b>Product Name:</b>	c-Fgr Monoclonal Antibody
<b>Catalog Number:</b>	N0150
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, IHC-p, IF (paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant fragment of human c-Fgr expressed in E. Coli.
<b>Concentration:</b>	
<b>Gene Name:</b>	FGR
<b>Protein Name:</b>	Tyrosine-protein kinase Fgr
<b>Other Name:</b>	FGR; SRC2; Tyrosine-protein kinase Fgr; Gardner-Rasheed feline sarcoma viral; v-fgr) oncogene homolog; Proto-oncogene c-Fgr; p55-Fgr; p58-Fgr; p58c-Fgr
<b>Human Swiss Prot No.:</b>	P09769
<b>Formulation:</b>	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
<b>Specificity:</b>	c-Fgr Monoclonal Antibody detects endogenous levels of c-Fgr protein.
<b>Storage and Stability:</b>	-20°C/1 year